If our standard configurations don't address your unique requirement, call us. We'll create a custom model just for you (subject to minimum quantities).




HC2TS27 HC5NMTS20





17"wx20"hx3.9"d



Dimensions	12.3"w x 26.5"h x 3.9"d	Din
Compartment Dimensions	Sanitizer Dispenser: 9.2"wx17.8"hx3.9"d Tissue/Glove compartment: 9"wx5.1"hx3.2"d	Co
Weight	12 lbs.	
Shipping Weight	17 lbs.	We
Description	(1) Sanitizer Station and single compartment for tissues or gloves.	
Options	Consumables not included with purchase of unit. Recommended materials: Purell ® TFX Touch-Free Dispenser (Model 2720-12) & Purell ® TFX 1200	Des
	ml Refill (Model 5456-04) Kimberly-Clark ® Professional Kleenex Facial Tissue 100ct. (Product Code 21400) Kimberly-Clark ® Safeskin Disposable Nitrile Gloves 100ct. Sizes Sm, Md, Lg & XL (Product Codes 55081, 55082, 55083 & 55084)	Ор
Mounting	Screw mount; template and all hardware for both drywall and masonry provided.	
Construction	Powder coated steel body with balance-laminated thermoplastic face, sides and top.	Мо
Finishes	See our standard material options on the product's web page: Or specify any custom finish.	Co
How to Specify	 Select model. Select face material. Determine quantity. 	Fin
Notes	*Prices subject to quantity	поч


Compartment Dimensions	Sanitizer Dispenser: 8"w x 18"h x 3.9"d Tissue/Glove compartment: 5.1"w x 9.6"h x 3.2"d Mask compartment: 3.7"w x 7"h x 3.2"d
Veight	13 lbs.
Shipping Weight	18 lbs.
Description	(1) Sanitizer Station, single compartment for tissues /gloves and single compartment for face masks.
Options	 Consumables not included with purchase of unit. Recommended materials: Purell ® TFX Touch-Free Dispenser (Model 2720-12) & Purell ® TFX 1200 ml Refill (Model 5456-04) Kimberly-Clark ® Professional Kleenex Facial Tissue 100ct. (Product Code 21400) Kimberly-Clark ® Safeskin Disposable Nitrile Gloves 100ct. Sizes Sm, Md, Lg & XL (Product Codes 55081, 55082, 55083 & 55084) Kimberly-Clark ® Procedure Mask 50ct. (Product Code 47080)
<i>l</i> ounting	Screw mount; template and all hardware for both drywall and masonry provided.
Construction	Powder coated steel body with balance-laminated thermoplastic face, sides and top.
inishes	See our standard material options on the product's web page: Or specify any custom finish.
How to Specify	 Select model. Select face material. Determine quantity.
lotes	*Prices subject to quantity

Dimensions	25.5"w x 21"h x 3.9"d
Compartment Dimensions	Sanitizer Dispenser: 8"w x 18"h x 3.9"d Tissue/Glove compartment: 5.1"w x 9.6"h x 3.2"d Mask compartment: 7"w x 3.7"h x 3.2"d Paper insert: 8.5"w x 11"h
Weight	21 lbs.
Shipping Weight	26 lbs.
Description	(1) Sanitizer Station, double compartment for tissues/gloves, single compartment for face masks and holder for paper insert.
Options	 Consumables not included with purchase of unit. Recommended materials: Purell ® TFX Touch-Free Dispenser (Model 2720-12) & Purell ® TFX 1200 ml Refill (Model 5456-04) Kimberly-Clark ® Professional Kleenex Facial Tissue 100ct. (Product Code 21400) Kimberly-Clark ® Safeskin Disposable Nitrile Gloves 100ct. Sizes Sm, Md, Lg & XL (Product Codes 55081, 55082, 55083 & 55084) Kimberly-Clark ® Procedure Mask 50ct. (Product Code 47080)
Mounting	Screw mount; template and all hardware for both drywall and masonry provided.
Construction	Powder coated steel body with balance-laminated thermoplastic face, sides and top.
Finishes	See our standard material options on the product's web page: Or specify any custom finish.
How to Specify	 Select model. Select face material. Determine quantity.
Notes	*Prices subject to quantity




HC4TGS21 HC5GTMS20 HC5GMG20





27.7"wx20"hx3.9"d




Dimensions	17"w×21"h×3.9"d	D
Compartment Dimensions	Sanitizer Dispenser: $8''w \times 18''h \times 3.9''d$ Tissue/Glove compartment: $5.1''w \times 9.6''h \times 3.2''d$	C D
Weight	14 lbs.	
Shipping Weight	19 lbs.	W
Description	(1) Sanitizer Station and double compartment for tissues/gloves.	S
Options	Consumables not included with purchase of unit. Recommended materials: • Purell ® TFX Touch-Free Dispenser (Model 2720-12) & Purell ® TFX 1200 ml Refill (Model 5456-04)	D
	 Kimberly-Clark ® Professional Kleenex Facial Tissue 100ct. (Product Code 21400) Kimberly-Clark ® Safeskin Disposable Nitrile Gloves 100ct. Sizes Sm, Md, Lg & XL (Product Codes 55081, 55082, 55083 & 55084) 	0
Mounting	Screw mount; template and all hardware for both drywall and masonry provided.	
Construction	Powder coated steel body with balance-laminated thermoplastic face, sides and top.	N
Finishes	See our standard material options on the product's web page: Or specify any custom finish.	С
How to Specify	 Select model. Select face material. Determine quantity. 	Fi
Notes	*Prices subject to quantity	Н

Compartment Dimensions	Sanitizer Dispenser: 8"w x 18"h x 3.9"d Tissue compartment: 5.1"w x 9.6"h x 3.2"d Mask compartment: 7"w x 3.7"h x 3.2"d Glove (x3) compartment: 5.1"w x 9.6"h x 3.2"d	C
Weight	21 lbs.	
Shipping Weight	26 lbs.	S
Description	(1) Sanitizer Station, single compartment for tissues, single compartment for face masks and triple compartment for glove boxes.	
Options	Consumables not included with purchase of unit. Recommended materials: • Purell ® TFX Touch-Free Dispenser (Model 2720-12) & Purell ® TFX 1200 ml Refill (Model 5456-04) • Kimberly-Clark ® Professional Kleenex Facial Tissue 100ct. (Product Code 21400) • Kimberly-Clark ® Safeskin Disposable Nitrile Gloves 100ct. Sizes Sm, Md, Lg & XL (Product Codes 55081, 55082, 55083 & 55084) • Kimberly-Clark ® Procedure Mask 50ct. (Product Code 47080)	N C
Mounting	Screw mount; template and all hardware for drywall and masonry wall provided	F
Construction	Powder coated steel body with balance-laminated thermoplastic face, sides and top.	Н
Finishes	See our standard material options on the product's web page: Or specify any custom finish.	
How to Specify	 Select model. Select face material. Determine quantity. 	N
Notes	*Prices subject to quantity	


Dimensions	25"wx20"hx3.9"d
Compartment Dimensions	Gown Compartment: $11.5''w \times 12.6''h \times 3.2''d$ Mask compartment: $7''w \times 3.7''h \times 3.2''d$ Glove (x3) Compartment: $9.6''w \times 5.1''h \times 3.2''d$
Weight	21 lbs.
Shipping Weight	26 lbs.
Description	Gown compartment, single compartment for masks and triple compartment for glove boxes.
Options	Consumables not included with purchase of unit. Recommended materials: Kimberly-Clark ® Safeskin Disposable Nitrile Gloves 100ct. Sizes Sm, Md, Lg & XL (Product Codes 55081, 55082, 55083 & 55084) Kimberly-Clark ® Procedure Mask 50ct. (Product Code 47080)
Mounting	Screw mount; template and all hardware for both drywall and masonry provided.
Construction	Powder coated steel body with balance-laminated thermoplastic face, sides and top.
Finishes	See our standard material options on the product's web page: Or specify any custom finish.
How to Specify	Select model. Select face material.
	3. Determine quantity.
Notes	*Prices subject to quantity *Prices subject to quantity



HC6GMGS20 HC4TG20








Dimensions	33.9"wx20"hx3.9"d
Compartment Dimensions	Sanitizer Dispenser: $8"w \times 18"h \times 3.9"d$ Gown compartment: $11.5"w \times 12.6"h \times 3.2"d$ Mask compartment: $7"w \times 3.7"h \times 3.2"d$ Glove (x3) compartment: $5.1"w \times 9.6"h \times 3.2"d$
Weight	25 lbs.
Shipping Weight	30 lbs.
Description	(1) Sanitizer Station, single compartment for masks, compartment for gowns and triple compartment for glove boxes.
Options	 Consumables not included with purchase of unit. Recommended materials: Purell ® TFX Touch-Free Dispenser (Model 2720-12) & Purell ® TFX 1200 ml Refill (Model 5456-04) Kimberly-Clark ® Professional Kleenex Facial Tissue 100ct. (Product Code 21400) Kimberly-Clark ® Safeskin Disposable Nitrile Gloves 100ct. Sizes Sm, Md, Lg & XL(Product Codes 55081, 55082, 55083 & 55084) Kimberly-Clark ® Procedure Mask 50ct. (Product Code 47080)
Mounting	Screw mount; template and all hardware for both drywall and masonry provided.
Construction	Powder coated steel body with balance-laminated thermoplastic face, sides and top.
Finishes	See our standard material options on the product's web page: Or specify any custom finish.
How to Specify	 Select model. Select face material. Determine quantity.
Notes	*Prices subject to quantity


Dimensions	10"wx20"hx3.9"d
Compartment Dimensions	Sanitizer Dispenser: 8"w x 18"h x 3.9"d
Weight	7 lbs.
Shipping Weight	12 lbs.
Description	(1) Sanitizer Station
Options	Consumables not included with purchase of unit. Recommended materials: • Purell ® TFX Touch-Free Dispenser (Model 2720-12) & Purell ® TFX 1200 ml Refill (Model 5456-04
Mounting	Screw mount; template and all hardware for both drywall and masonry provided.
Construction	Powder coated steel body with balance-laminated thermoplastic face, sides and top.
Finishes	See our standard material options on the product's web page: Or specify any custom finish.
How to Specify	 Select model. Select face material. Determine quantity.
Notes	*Prices subject to quantity

Dimensions	18.4"w×21"h×3.9"d
Compartment Dimensions	Tissue compartment: 5.1"w x 9.6"h x 3.2"d Glove (x3) compartment: 5.1"w x 9.6"h x 3.2"d
Weight	18 lbs.
Shipping Weight	23 lbs.
Description	(2) Single compartments for tissues and triple compartment for glove boxes.
Options	 Consumables not included with purchase of unit. Recommended materials: Kimberly-Clark ® Professional Kleenex Facial Tissue 100ct. (Product Code 21400) Kimberly-Clark ® Safeskin Disposable Nitrile Gloves 100ct. Sizes Sm, Md, Lg & XL (Product Codes 55081, 55082, 55083 & 55084)
Mounting	Screw mount; template and all hardware for both drywall and masonry provided.
Construction	Powder coated steel body with balance-laminated thermoplastic face, sides and top.
Finishes	See our standard material options on the product's web page: Or specify any custom finish.
How to Specify	 Select model. Select face material. Determine quantity.
Notes	*Prices subject to quantity



HC2G17 HC4G12 HC1W3







Compartment Dimensions Clower (x3) compartment 5.1°w x 9.8°h x 3.2°d Compartment Dimensions Glove (x3) compartment 5.1°w x 9.8°h x 3.2°d Weight 7 lbs. Weight 8 lbs. Shipping Weight 12 lbs. Shipping Weight 13 lbs. Description Triple compartment glove holder. Description (3) Single glove compartments Options Consumables not included with purchase of unit. Recommended materials: Included with	Dimensions	11.5"w×16.7"h×3.9"d	Dimensions	17.2"w×11.2"h×3.9"d
Shipping Weight 12 lbs. Shipping Weight 13 lbs. Description Triple compartment glove holder. Description (3) Single glove compartments Consumables not included with purchase of unit. Recommended materials: Kimberly-Clark (® Safeskin Disposable Nitrile Gloves 100ct. Sizes Sm, Md, Lg & XL (Product Codes 55081, 55082, 55083 & 55084) Mounting Screw mount; template and all hardware for both drywall and masonry provided. Construction Powder coated steel body with balance-laminated thermoplastic face, sides and top. Finishes See our standard material options on the product's web page: Or specify any custom finish. How to Specify 1. Select model. 2. Select face material. 3. Determine quantity. Shipping Weight 13 lbs. Description (3) Single glove compartments Consumables not included with purchase of unit. Recommended materials: Consumables not included with purchase of unit. Recommended materials: Consumables not included with purchase of unit. Recommended materials: Consumables not included with purchase of unit. Recommended materials: Consumables not included with purchase of unit. Recommended materials: Consumables not included with purchase of unit. Recommended materials: Consumables not included with purchase of unit. Recommended materials: Consumables not included with purchase of unit. Recommended materials: Consumables not included with purchase of unit. Recommended materials: Consumables not included with purchase of unit. Recommended materials: Consumables not included with purchase of unit. Recommended materials: Consumables not included with purchase of unit. Recommended insertials: Consumables not included with purchase of unit. Recommended insertials: Consumables not included with purchase of unit. Recommended insertials: Consumables not included with purchase of unit. Recommended insertials: Consumables not included with purchase of unit. Recommended insertials: Consumables not included with purchase of unit. Recommended insertials: Consumables not included with purchase of u	•	Glove (x3) compartment: 5.1"w x 9.6"h x 3.2"d	•	Glove (x3) compartment: 5.1"w x 9.6"h x 3.2"d
Description Triple compartment glove holder. Consumables not included with purchase of unit. Recommended materials:	Weight	7 lbs.	Weight	8 lbs.
Consumables not included with purchase of unit. Recommended materials: Vimberly-Clark ® Safeskin Disposable Nitrile Gloves 100ct. Sizes Sm, Md, Lg & XL (Product Codes 55081, 55082, 55083 & 55084) Mounting Screw mount; template and all hardware for both drywall and masonry provided. Mounting Powder coated steel body with balance-laminated thermoplastic face, sides and top. Finishes See our standard material options on the product's web page: Or specify any custom finish. How to Specify 1. Select model. 2. Select face material. 3. Determine quantity. Consumables not included with purchase of unit. Recommended materials: Consumables not included with purchase of unit. Recommended materials: Vimberly-Clark ® Safeskin Disposable Nitrile Gloves 100ct. Sizes Sm, Md, Lg & XL (Product Codes 55081, 55082, 55083 & 55084) Mounting Screw mount; template and all hardware for drywall and masonry wall provided Construction Powder coated steel body with balance-laminated thermoplastic face, sides and top. Finishes See our standard material options on the product's web page: Or specify any custom finish. How to Specify 1. Select model. 2. Select face material. 3. Determine quantity.	Shipping Weight	12 lbs.	Shipping Weight	13 lbs.
 Options Kimberly-Clark ® Safeskin Disposable Nitrile Gloves 100ct. Sizes Sm, Md, Lg & XL (Product Codes 55081, 55082, 55083 & 55084) Options Kimberly-Clark ® Safeskin Disposable Nitrile Gloves 100ct. Sizes Sm, Md, Lg & XL (Product Codes 55081, 55082, 55083 & 55084) Mounting Screw mount; template and all hardware for both drywall and masonry provided. Construction Powder coated steel body with balance-laminated thermoplastic face, sides and top. Finishes See our standard material options on the product's web page: Or specify any custom finish. How to Specify 1. Select model. 2. Select face material. 3. Determine quantity. 	Description	Triple compartment glove holder.	Description	(3) Single glove compartments
Powder coated steel body with balance-laminated thermoplastic face, sides and top. Finishes See our standard material options on the product's web page: Or specify any custom finish. Finishes 1. Select model. 2. Select face material. 3. Determine quantity. Powder coated steel body with balance-laminated thermoplastic face, sides and top. Finishes See our standard material options on the product's web page: Or specify any custom finish. Finishes See our standard material options on the product's web page: Or specify any custom finish. How to Specify 1. Select model. 2. Select face material. 3. Determine quantity.	Options	· Kimberly-Clark ® Safeskin Disposable Nitrile Gloves 100ct. Sizes Sm, Md,	Options	· Kimberly-Clark ® Safeskin Disposable Nitrile Gloves 100ct. Sizes Sm, Md,
and top. Finishes See our standard material options on the product's web page: Or specify any custom finish. Finishes See our standard material options on the product's web page: Or specify any custom finish. Finishes See our standard material options on the product's web page: Or specify any custom finish. How to Specify 1. Select model. 2. Select face material. 3. Determine quantity. 1. Select model. 2. Select face material. 3. Determine quantity.	Mounting		Mounting	
custom finish. How to Specify 1. Select model. 2. Select face material. 3. Determine quantity. How to Specify 2. Select face material. 3. Determine quantity. custom finish. How to Specify 2. Select model. 2. Select face material. 3. Determine quantity.	Construction	, , ,	Construction	. , ,
2. Select face material. 2. Select face material. 3. Determine quantity. 2. Select face material. 3. Determine quantity.	Finishes		Finishes	
Notes *Prices subject to quantity *Prices subject to quantity	How to Specify	2. Select face material.	How to Specify	2. Select face material.
	Notes	*Prices subject to quantity	Notes	*Prices subject to quantity

Dimensions	6"w×11.5"h×3.94"d
Compartment Dimensions	N/A
Weight	1 lbs.
Shipping Weight	4 lbs.
Description	Small, Disinfectant wipe holder.
Options	Consumables not included with purchase of unit.
Mounting	Screw mount; template and all hardware for both drywall and masonry provided.
Construction	Powder coated steel bracket on balance-laminated wood composite core backer.
Finishes	See our standard material options on the product's web page: Or specify any custom finish.
How to Specify	 Select model. Select backer material. Determine quantity.
Notes	*Prices subject to quantity



HC1W4 HC1W5





Dimensions	7"w×11.5"h×5.25"d	Dimensions	8"wx11.5"hx6.25"d
Compartment Dimensions	N/A	Compartment Dimensions	N/A
Weight	1 lbs.	Weight	1 lbs.
Shipping Weight	4 lbs.	Shipping Weight	4 lbs.
Description	Medium, Disinfectant wipe holder.	Description	Large, Disinfectant wipe holder.
Options	Consumables not included with purchase of unit.	Options	Consumables not included with purchase of unit.
Mounting	Screw mount; template and all hardware for both drywall and masonry provided.	Mounting	Screw mount; template and all hardware for both drywall and masonry provided.
Construction	Powder coated steel bracket on balance-laminated wood composite core backer.	Construction	Powder coated steel bracket on balance-laminated wood composite core backer.
Finishes	See our standard material options on the product's web page: Or specify any custom finish.	Finishes	See our standard material options on the product's web page: Or specify any custom finish.
How to Specify	 Select model. Select backer material. Determine quantity. 	How to Specify	 Select model. Select backer material. Determine quantity.
Notes	*Prices subject to quantity	Notes	*Prices subject to quantity